


**Bondable Traction Hook**

Traction Hooks encourage and direct the eruption of severely malpositioned or impacted teeth that have had a portion of the crown surgically exposed. Eyelets are securely soldered to 3.56 mm (.140") diameter laminated mesh bonding bases. Base can be trimmed or recontoured, if necessary, to fit whatever portion of the tooth that is exposed.

Square traction hook has lower profile than round hook, which increases patient comfort, yet provides 30% more room to facilitate threading of steel or elastomeric ligatures.

Round traction hook with prefabricated ligating chain can be used when tooth being bonded is deeply imbedded in gingival tissue. The 34 mm chain is formed from 0.31 mm (.012") ligature wire and can be easily cut to required length at each appointment as tooth continues to erupt.




224-011 Round Traction Hook  
224-012 Round Traction Hook with Chain  
224-013 Square Traction Hook  
Package of 10

Heat Sterilize

**Reverse Palatal Arch Bar**


Palatal Arch Bars provide excellent control of upper first molars that require stabilization, rotation, forcing, or expansion. Reverse (mesially oriented) adjustment loop is less likely to cause a gag reflex.

Molars can be rotated to the buccal by compressing the palatal loop during insertion, or molar stability can be maintained by inserting a passive Palatal Arch Bar. Bars are designed to fit 0.91 mm (.036") lingual sheaths and are removable. Ligating Palatal Bars to the lingual sheath is recommended.



mm	Order No.	mm	Order No.
37	352-037	49	352-049
40	352-040	52	352-052
43	352-043	Kit*	352-100
46	352-046		

\*Kit contains 30 arches of assorted sizes and a measuring device.  
Individual Sizes: Package of 5



Heat Sterilize

Molar sheaths are available on page 83.